


Glidewire Advantage™

Guide Wire

OFFERS YOU AN **ADVANTAGE** FROM ACCESS TO TREATMENT



High durability¹

Shape retention capability of nitinol core provides high durability.

	GWA 0.014"	GWA 0.018"
Before		
After		


Test Protocol This test evaluated the distal tip of the guidewire after turning and bending it.



Nitinol core

Great crossability¹


Optimal lubricity design in the first 25cm distal portion of the wire leads to great crossability.



Test Protocol Measure the sliding resistance of hydrophilic coating of each guidewire

Improved steerability¹


The spiral PTFE structure acts as a grip for better handling and navigational control.



Outside shaft grip

Stable device support¹

A stiff proximal core provides device support even in complex anatomies.



Stiff proximal nitinol core